1. Intended Use

For *In Vitro* Diagnostic Use

The IDS-iSYS 1,25 VitD^{Xp} assay is an *in vitro* diagnostic device intended for the quantitative determination of 1,25-dihydroxyvitamin D [1,25(OH)₂D] in human serum on the IDS-iSYS Multi-Discipline Automated System. Results are to be used as an aid in the assessment of vitamin D sufficiency.

Rx only.

2. Summary and Explanation

There are two forms of vitamin D: vitamin D₃ and vitamin D₂. Vitamin D₃, cholecalciferol, is the naturally occurring vitamin D that is produced in the skin after the exposure of 7-dehydrocholesterol to solar ultraviolet radiation¹. Vitamin D₂ is manufactured through the ultraviolet irradiation of ergosterol from yeast. Both are used in vitamin D supplements. The vitamin D compound is biologically inactive but enters the circulation and is hydroxylated in the liver to 25-hydroxyvitamin D [25(OH)D], which is used to determine a patient's vitamin D status.

In the kidney, 25(OH)D is further hydroxylated to produce the biologically active metabolite, 1,25(OH)₂D¹. 1,25(OH)₂D is one of the major regulators of calcium and phosphate metabolism, stimulating intestinal calcium absorption and increasing bone resorption. 1,25(OH)₂D also inhibits parathyroid hormone (PTH) production both by direct action on the parathyroid glands and indirectly by raising serum calcium levels. 1,25(OH)₂D production is itself stimulated by parathyroid hormone (PTH), thus providing an effective control loop¹. Fibroblast growth factor 23 [FGF-23], secreted from the bone, causes the sodium-phosphate co-transporter to be internalised by the cells of the kidney and small intestine and also suppresses 1,25(OH)₂D synthesis².

In secondary hyperparathyroidism, a disease outside of the parathyroids causes the parathyroid glands to become enlarged and hyperactive. It is usually caused by kidney failure, a problem where the kidney is unable to cleanse the blood of phosphorus produced by the body and unable to make enough vitamin D, specifically 1,25(OH)₂D - the active form of vitamin D. The build-up of phosphorus leads to low levels of calcium in the blood, which stimulates the parathyroid glands to increase parathyroid hormone (PTH) production leading to enlargement of the glands. As the disease progresses, the parathyroid glands no longer respond normally to calcium and Vitamin D. The rationale for direct activated vitamin D therapy in CKD is to slow the progression of secondary hyperparathyroidism. Clinical practice guidelines such as the 'Kidney Disease Outcomes Quality Initiatives (KDOQI)' and 'Kidney Disease: Improving Global Outcomes (KDIGO)' recommend activated vitamin D therapeutic regimens for CKD patients^{3,4}.

3. Method Description

The IDS-iSYS 1,25VitD^{Xp} kit is a complete test system involving immunoextraction of 1,25-dihydroxyvitamin D in human serum and quantitative determination of 1,25(OH)₂D levels on the IDS-iSYS Multi-Discipline Automated System (System).

In the cuvette 1, 218 µL of patient sample is delipidated then incubated with a specific anti-1,25(OH)₂D antibody coated magnetic particles. After incubation, the magnetic particles are washed and the 1,25(OH)₂D in the patient sample is eluted. 180 µL of the extracted 1,25(OH)₂D sample is transferred to cuvette 2 where the immunoassay, based on chemiluminescence technology, takes place.

In the cuvette 2, 180 µL purified samples are incubated with biotinylated sheep anti-1,25(OH)₂D antibody. Following the first incubation step, an acridinium labelled 1,25(OH)₂D conjugate is added followed by a subsequent incubation step. Streptavidin labelled magnetic particles are then added prior to a third incubation step. The magnetic particles are captured using a magnet and a wash step performed to remove any unbound 1,25(OH)₂D. Trigger reagents are added; the resulting light emitted by the acridinium label is inversely proportional to the concentration of 1,25(OH)₂D in the original sample.

4. Warnings and Precautions

The IDS-iSYS 1,25VitD^{Xp} is for *in vitro* diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in these Instructions For Use (IFU). Immunodiagnostic Systems Limited (IDS) will not be held responsible for any loss or damage (except as required by statute), howsoever caused, arising out of non-compliance with the instructions provided.

CAUTION: This kit contains material of animal origin. Handle kit reagents as if capable of transmitting an infectious agent. Appropriate precautions and good laboratory practice must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Human materials

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV I and II), Hepatitis B surface antigen, antibody to Hepatitis C, and found to be negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled according to Biosafety Level 2.

Reagents containing Sodium Azide

Some reagents in this kit contain sodium azide (NaN₃) <0.1 % (w/w) which may react with lead, copper or brass plumbing to form highly explosive metal azides. On disposal, flush with large volumes of water to prevent azide build up.

Classification under CLP:

Acute Tox 4
STOT SE 1
Eye Irrit. 2
Skin Irrit 2

Hazard statements:

H302+312: Harmful if swallowed or in contact with skin.
H315: Causes skin irritation.
H319: Causes serious eye irritation.
H370: Causes damage to organs .

Precautionary statements:

P260: Do not breathe vapours.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P301+312: IF SWALLOWED: Call a POISON CENTER/doctor/ if you feel unwell.
P302+350: IF ON SKIN: Gently wash with plenty of soap and water.
P302+352: IF ON SKIN: Wash with plenty of water/.
P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

5. Handling Precautions

Apart from the calibrators, which are lyophilised, the reagents provided in the kit are ready to use. Refer to the Calibrator section of the procedure for reconstitution methodology.

Before a new cartridge is loaded onto the System, the magnetic particle container requires mixing by the operator with a brisk rotation motion. This will re-suspend the magnetic particles that have settled during shipment. Care must be taken to avoid foam formation in the reagent cartridge.

Shelf Life and Storage of Reagents

Store the cartridge and calibrators in an **upright** position in the dark at 2 - 8 °C. **Do not freeze** the cartridge or the calibrators.

Reagent shelf life	Immunextraction Cartridge	Immunoassay Cartridge	Calibrators
Before opening at 2 - 8 °C	To the expiry date		
Cartridge, after opening at 2 - 8 °C	N/A	56 Days	N/A
Calibrators, after reconstitution at 2 - 8 °C	N/A	N/A	6 Hours
On board the System *	28 Days	49 Days	3 Hours

* Continuous on board stability

6. Sample Collection and Storage

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) samples.

- Samples should be separated as soon as possible after collection.
- Following separation, samples can be stored for a maximum of 24 hours at room temperature (18 - 25°C) or up to 7 days at 2-8°C. If longer time is required between separation and testing, the samples should be stored at -20°C.
- Samples can be stored at -20°C for up to 2.5 months, following which they should be stored at -80°C.
- Samples may undergo up to 2 freeze thaw cycles and should be completely thawed and mixed well prior to loading onto the instrument.
- Samples are stable on-board the analyser for a maximum of 3 hours.

Note:

- Samples containing particulate matter must be centrifuged before performing the assay. Do not use heat-inactivated samples.
- To minimise possible evaporation effects, samples, calibrators and controls should be measured within 3 hours of being placed on the System.
- Each laboratory should follow the guidelines or requirements of local, state, and/or federal regulations or accrediting organizations to establish its own specimens handling and storage stability. For guidance on appropriate practices, please refer to the CLSI GP44-A4, "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests".
- Before performing assays, make sure that samples, calibrators and controls are at room temperature (18 - 25°C).
- Some commercially available sample collection tubes may affect the results of testing in particular cases.
- It is recommended to follow the instructions of the tube manufacturer especially when processing samples in primary tubes.
- Centrifuged samples with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified samples without the lipemic material.
- **The samples MUST equilibrate to room temperature (18 - 25°C) before loading onto the System in order to optimise the delipidation step.**

7. Materials

Materials Provided

Immunoextraction Cartridge, 2 cartridges of 50 tests

MPA1	Magnetic particles coated with monoclonal antibody to 1,25(OH) ₂ D, 1 vial, 2.8 mL.
DEL	Delipidation reagent, a solution of dextran sulphate and magnesium chloride, 1 vial, 2.6 mL.
DIS	Displacer reagent, a proprietary reagent for dissociating 1,25(OH) ₂ D from binding proteins with sodium azide as preservative (<0.1%), 1 vial, 18.7 mL.
ELU	Elution reagent, a solution of sodium hydroxide, 1 vial, 13.0 mL.
NET	Neutralisation reagent, a solution of citric acid, 1 vial, 3.4 mL.
BUF	Assay buffer, a MOPS buffer containing bovine serum albumin, 1vial, 9 mL per vial.

Immunoassay Cartridge, 1 cartridge of 100 tests

MPS3	Magnetic particles coated with streptavidin in a phosphate buffer containing bovine serum albumin with sodium azide as a preservative (<0.1%), 1 vial, 2.6 mL.
CONJ	Conjugate, 1,25(OH) ₂ D labelled with an acridinium ester derivative, in a phosphate buffer containing bovine serum albumin with sodium azide as preservative (<0.1%), 1 vial, 9.6 mL.
Ab-BIOT	Antibody-Biotin, anti-1,25(OH) ₂ D polyclonal antibody labelled with biotin in a phosphate buffer containing sheep proteins and sodium azide as preservative (<0.1%), 1 vial, 13 mL.
BUFD	Wash buffer, 1% proprietary detergent in PBS with sodium azide as a preservative (<0.1%), 1 vial, 35 mL per vial.

Calibrators

CAL A	Lyophilised MOPS buffer containing bovine serum albumin, 1,25(OH) ₂ D with sodium azide as preservative (<1.0%), 3 each of 2 concentration levels, 1.2 mL per bottle.
CAL B	

Mini CD

Contains IFU for IDS-iSYS reagents, CRY files and Certificate of Analysis.

Materials Required But Not Provided

System

IS-310400	IDS-iSYS Multi-Discipline Automated System
IS-2030	IDS-iSYS 1,25 VitD ^{xp} Control Set
IS-CC100	IDS-iSYS Cuvettes
IS-CS100	IDS-iSYS System Liquid (Syst I)
IS-CW100	IDS-iSYS Wash Solution (Wash S)
IS-CT100	IDS-iSYS Trigger Set
IS-6010	IDS-iSYS Cartridge Check System (CCS)
IS-CSC105	Sample Cups (500 µL)
IS-DW225	Disposable waste bags
IS-DS100	D-SORB solution

Optional equipment / materials

IS-2035	IDS-iSYS 1,25 VitD ^{xp} Calibration Verifiers
IS-10DA	IDS-iSYS Diluent A
IS-XP01	IDS-iSYS XPrep
	Disposable polypropylene 2 mL, 10.8 mm diameter conical skirted base, screw cap micro tubes and screw cap with O-ring (Sarstedt 72.609 and 65.716 or equivalent)
	Precision Pipetting Devices

8. Assay Procedure

Immunoextraction Reagent Cartridge

The reagents provided in the cartridge are ready to use. Before a new cartridge is loaded on board the System, the magnetic particles container must be completely re-suspended. The magnetic particles can be prepared following using one of following options, avoid foam formation:

Manually:

- Mix the magnetic particles vial briskly between the palms of the hands in a back and forth motion until the colour of the suspended particles has changed to brown.
- Gently and carefully invert the magnetic particles vial until all settled magnetic particles in the bottom of the vial have re-suspended.
- Visually check the bottom of the magnetic particles vial to confirm that all settled magnetic particles have completely dispersed.
- Repeat the inversion as necessary until the magnetic particles are completely re-suspended.
- Ensure that there is no foam formation in the cartridge reagents.

End-over end test/blood tube rotator:

- Gently and carefully remove the particle vial from the plastic clip.
- Mount the particles vial on the test/blood holder.
- Rotate the vial at 10-20 revolutions per minute for 15 minutes at room temperature (18-25°C).
- Visually check the bottom of the magnetic particles vial to confirm that all settled magnetic particles have completely dispersed.
- Repeat the inversion as necessary until the magnetic particles are completely re-suspended.
- Ensure that there is no foam formation in the cartridge reagents.

IDS-iSYS XPrep (IS-XP01):

- Refer to the IDS-iSYS XPrep user manual for usage instructions.

The barcode is read when the cartridge is loaded on the reagent tray. If the label cannot be read by the System's barcode reader, a manual procedure exists to enter the barcode data (see IDS-iSYS User Manual). Once the cartridge is on the reagent tray, the System automatically performs the mixing of magnetic particles to maintain homogeneity.

Do not perform the cartridge level scan on the IDS-iSYS 1,25 VitD^{xp} Immunoextraction cartridge as this will result in less than 50 tests being obtained from the cartridge.

The cartridge should be loaded on the reagent tray for at least 60 minutes before starting the assay. Leave the immunoextraction reagent cartridge on-board the System until the cartridge has been fully depleted.

Immunoassay Reagent Cartridge

The reagents provided in the cartridge are ready to use. Before a new cartridge is loaded on board the System, mix the magnetic particles container with a brisk rotation motion between the palms of the hands. Avoid foam formation.

IDS-iSYS XPrep (IS-XP01):

- Refer to the IDS-iSYS XPrep user manual for usage instructions

The barcode is read when the cartridge is loaded on the reagent tray. If the label cannot be read by the System's barcode reader, a manual procedure exists to enter the barcode data (see IDS-iSYS User Manual). Once the cartridge is on the reagent tray, the System automatically performs the mixing of magnetic particles to maintain homogeneity.

Do not perform the cartridge level scan on the IDS-iSYS 1,25 VitD^{xp} Immunoassay cartridge as this will result in less than 100 tests being obtained from the cartridge.

The cartridge should be loaded on the reagent tray for at least 40 minutes before starting the assay. If the cartridge is removed from the reagent tray, store the cartridge upright at 2 - 8°C in the dark.

Calibrators

The IDS-iSYS 1,25 VitD^{xp} assay has been standardised against in-house reference standards which has been aligned to the validated 1,25-dihydroxyvitamin D LC-MS/MS.

The IDS-iSYS 1,25 VitD^{xp} calibrators are lyophilised. Prepare the calibrators following below procedure:

- Ensure the calibrators equilibrate to room temperature (18 - 25°C).
- Immediately before use, add 1.2 mL of distilled or deionised water to each bottle.
- Replace the stopper. Allow the vials to stand for 10 minutes at room temperature (18-25°C).
- Vortex for 2-3 seconds; invert gently to ensure complete reconstitution.
- Transfer approximately 750 µL for triplicate calibration (180 µL per replicate) of calibrators to polypropylene 2 mL, conical skirted base, screw cap micro tubes and place onto the System within 15 minutes of reconstitution and measure within 3 hours. The calibrators do not go through the purification procedure
- Avoid the formation of foam during each step.
- Proceed according to the instructions in the IDS-iSYS User Manual.

Calibration Frequency

A new calibration is required:

- Upon loading each new lot of cartridges.
- When the control values do not fall within the defined ranges.
- When the calibration interval of 21 days has expired.
- After System service.

Verification of the calibration is automatic and managed by the System.

System Calibration

The two IDS-iSYS 1,25 VitD^{xp} calibrators are required to perform the adjustment of the master curve. The calibrators are supplied with the kit and calibrators from another lot must not be used. All levels of IDS-iSYS 1,25 VitD^{xp} Control Set (IS-2030) MUST be measured at the same time as the calibrators to validate the master curve adjustment.

All data required for the calibration of the cartridge can be found on the mini CD. Use calibrator levels A and B to adjust the master curve to the reagents on board the System. Check for the presence of an IDS-iSYS 1,25 VitD^{xp} Immunoextraction and an Immunoassay cartridge on the reagent tray and the availability of the cartridge master curve in the database. If the data for the lot of calibrators is not available on board the System, load the data using the mini CD provided with the reagent kit.

Start the immunoassay calibration on the System according to the System User Manual. The calibration is carried out using 3 replicates of each level. One replicate may be removed to meet the calibration requirements. All control levels must also be measured in duplicate to calibrate the assay. Verify and approve the calibration according to the calibration status displayed in the calibration window and discard the remaining volume of calibrators and controls from the sample tray after use.

9. Quality Control

Use the IDS-iSYS 1,25 VitD^{xp} Control Set for quality control. To ensure validity of results at least two controls with varying levels of 1,25(OH)₂D should be measured. Other suitable control material can be used in addition to the IDS-iSYS 1,25 VitD^{xp} Control Set. Controls should be tested at (or near) the beginning of every run containing patient samples and also during every calibration in accordance with local, state and/or federal regulations or accreditation requirements and your laboratory's quality procedure. Refer to the IDS-iSYS 1,25 VitD^{xp} Control Set IFU for preparation and handling instructions.

Determination of Sample 1,25-dihydroxyvitamin D (1,25(OH)₂D) levels

Process samples according to the IDS-iSYS System User Manual.

10. Calculation of Results

The concentration of 1,25(OH)₂D within each sample is calculated automatically. Display of the concentrations (screen or printed) is produced upon user request.

The IDS-iSYS 1,25 VitD^{xp} assay uses a 4-parameter logistic curve fit (4PL) to calculate the 1,25(OH)₂D concentrations.

To convert results to SI units: pmol/L = pg/mL x 2.4

11. Measurement Range (Reportable Range)

The reportable range of the assay is 7.5 - 150 pg/mL (18 – 360 pmol/L). Any value that reads below 7.5 pg/mL (18 pmol/L) should be reported as “< 7.5 pg/mL” (“< 18 pmol/L”).

12. Dilution

Samples with 1,25(OH)₂D concentrations above the reportable range should be diluted manually with IDS-iSYS Diluent A (IS-10DA) in a ratio of 1 in 2. The results for diluted samples must be multiplied by the 2x dilution factor.

13. Limitations of Use

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient’s clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁵. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The following substances do not interfere in the IDS-iSYS 1,25 VitD^{xp} assay when the concentrations presented in the following table are below the stated threshold.

Potentially Interfering Agent	Threshold Concentration
Albumin	7.5 g/dL
Bilirubin, conjugated	40 mg/dL
Bilirubin, unconjugated	40 mg/dL
Biotin	0.1 mg/dL
Cholesterol, total	300 mg/dL
HAMA	1000 ng/mL
Haemoglobin	500 mg/dL
Rheumatoid Factor	7000 IU/mL
Triglycerides	3000 mg/dL
Acetaminophen	9 mg/dL
Ibuprofen	50 mg/dL
Carbamazepine	127 µmol/L
Phenytoin	198 µmol/L
Acetylsalicylic Acid	65 mg/dL
Salicylic Acid	60 mg/dL

14. Expected Values

The following ranges were determined using the IDS-iSYS 1,25 VitD^{xp} assay and are provided for information only.

Samples from 121 apparently healthy male and female adults aged 19 – 69 years of age from mix ethnic backgrounds (34% dark-skinned and 66% light-skinned) were assayed in the IDS-iSYS 1,25 VitD^{xp}. Samples were collected during the winter (36%) and summer (64%) from individuals with normal intact PTH, calcium, phosphate, and TSH values from the northern and southern regions of the United States. The 95 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

Obtained normal adults range: 15.2 to 90.1 pg/mL (36.5 to 216.2 pmol/L)

15. Performance Data

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1 Sensitivity

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation”.

Sensitivity	1,25(OH) ₂ D Concentration	
	pg/mL	pmol/L
LoB (Limit of Blank)	2.5	6.0
LoD (Limit of Detection)	4.2	10.1
LoQ (Limit of Quantitation)	7.5	18.0

15.2 Precision

Precision was evaluated in accordance with a modified protocol based on CLSI EP-5A2, "Evaluation of Precision Performance of Quantitative Measurement Methods". A total of 6 samples were assayed in duplicate, twice per day for 7 days using one lot of reagents on one system.

Sample	n	Mean		Repeatability SD		CV%	Total SD		CV%
		pg/mL	pmol/L	pg/mL	pmol/L		pg/mL	pmol/L	
Serum 1	28	15.0	36.0	3.0	7.2	10.8	3.1	7.4	11.7
Serum 2	28	27.7	66.5	1.8	4.3	4.8	2.3	5.5	9.8
Serum 3	28	42.9	103.0	4.6	11.0	9.8	10.4	25.0	10.2
Serum 4	28	87.3	209.5	3.2	7.7	8.4	8.8	21.1	13.0
Control 1	28	24.2	58.1	7.4	17.8	7.6	11.4	27.4	9.5
Control 2	28	80.4	193.0	11.8	28.3	4.0	19.9	47.8	10.9

15.3 Linearity

Linearity was evaluated based on CLSI EP-6A, "Evaluation of the Linearity of Quantitative Measurement Procedures". Samples were prepared by diluting a high patient sample with a low patient sample prior to assay. The linear regression of the Observed concentrations versus the Expected concentrations (Expected = Slope x Observed ± Intercept) is:

Slope	Intercept		Regression coefficient (R ²)
	pg/mL	pmol/L	
1.08	-6.5	-15.6	0.99

15.4 Method Comparison

The IDS-iSYS 1,25 VitD^{xp} assay was compared against the commercially available quantitative IDS-iSYS 1,25-Dihydroxy Vitamin D assay, following CLSI EP-9A2, "Method Comparison and Bias Estimation Using Patient Samples". A total of 105 samples, selected to represent a wide range of 1,25(OH)₂D concentrations [11.4 - 150 pg/mL (27.4 – 360 pmol/L)], were assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope	Intercept		Correlation coefficient (r)
		pg/mL	pmol/L	
105	0.97	6.0	14.4	0.94

15.5 Specificity

The following exogenous synthetic vitamin D metabolites or vitamin D drug compounds were spiked into serum samples and measured with the IDS-iSYS 1,25 VitD^{xp} assay.

Analyte	Cross-Reactivity (%)
1,25(OH) ₂ D ₃	97%
1,25(OH) ₂ D ₂	72%
1,24,25(OH) ₃ D ₃	106.7%
25-OH D ₃	0.04%
25-OH D ₂	0.16%
C3-epimer 25(OH)D ₃	0.00%
24,25(OH) ₂ D ₃	0.06%
24,25(OH) ₂ D ₂	0.00%
25,26(OH) ₂ D ₃	0.52%
Alfacalcidol	0.01%
Cholecalciferol (Vitamin D3)	0.00%
Ergocalciferol (Vitamin D2)	0.00%
Sensipar	0.00%
Zemplar (paracalcitol)	5.70%

16. Symbols used

Catalogue Number

In Vitro Diagnostic Device

Manufacturer

Applied in accordance with directive 98/79/EC

Rx Only

Caution: US Federal Law restricts this device to sale by or on the order of a (licensed healthcare practitioner)

17. Bibliography

1. Holick MF. Vitamin D deficiency. *N Engl J Med* 2007;357:266-281.
2. Shimada T et al. FGF-23 is a potent regulator of vitamin D metabolism and phosphate homeostasis. *J Bone Miner Res* 2004;19:429-35.
3. K/DOQI clinical practice guidelines for bone metabolism and disease in chronic kidney disease. *Am J Kidney Dis* 2003; 42:Suppl 3:S1-S201.
4. Kidney Disease: Improving Global Outcomes (KDIGO) CKD-MBD Work Group. KDIGO clinical practice guideline for the diagnosis, evaluation, prevention, and treatment of chronic kidney disease-mineral and bone disorder (CKD-MBD). *Kidney Int.* 2009;76 (Suppl 113): S1-S130.
5. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idsplc.com • www.idsplc.com

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Hazard statements:

Precautionary statements:

H302+312: Harmful if swallowed or in contact with skin.	P260: Do not breathe vapours.
H315: Causes skin irritation.	P280: Wear protective gloves/protective clothing/eye protection/face protection.
H319: Causes serious eye irritation.	P301+312: IF SWALLOWED: Call a POISON CENTER/doctor/ if you feel unwell.
H370: Causes damage to organs .	P302+350: IF ON SKIN: Gently wash with plenty of soap and water.
	P302+352: IF ON SKIN: Wash with plenty of water/.
	P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

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Skin Irrit 2

Hazard statements:

H302+312: Harmful if swallowed or in contact with skin.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H370: Causes damage to organs .

Precautionary statements:

P260: Do not breathe vapours.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P301+312: IF SWALLOWED: Call a POISON CENTER/doctor/ if you feel unwell.

P302+350: IF ON SKIN: Gently wash with plenty of soap and water.

P302+352: IF ON SKIN: Wash with plenty of water/.

P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

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1. Holick MF. Vitamin D deficiency. *N Engl J Med* 2007;357:266-281.
2. Shimada T et al. FGF-23 is a potent regulator of vitamin D metabolism and phosphate homeostasis. *J Bone Miner Res* 2004;19:429-35.
3. K/DOQI clinical practice guidelines for bone metabolism and disease in chronic kidney disease. *Am J Kidney Dis* 2003; 42:Suppl 3:S1-S201.
4. Kidney Disease: Improving Global Outcomes (KDIGO) CKD-MBD Work Group. KDIGO clinical practice guideline for the diagnosis, evaluation, prevention, and treatment of chronic kidney disease-mineral and bone disorder (CKD-MBD). *Kidney Int.* 2009;76 (Suppl 113): S1-S130.
5. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

 Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Τηλ.: +44 191 519 0660 • Φαξ: +44 191 519 0760
Email: info.uk@idsplc.com • www.idsplc.com

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17. Bibliográfia

1. Holick MF. Vitamin D deficiency. *N Engl J Med* 2007;357:266-281.
2. Shimada T et al. FGF-23 is a potent regulator of vitamin D metabolism and phosphate homeostasis. *J Bone Miner Res* 2004;19:429-35.
3. K/DOQI clinical practice guidelines for bone metabolism and disease in chronic kidney disease. *Am J Kidney Dis* 2003; 42:Suppl 3:S1-S201.
4. Kidney Disease: Improving Global Outcomes (KDIGO) CKD-MBD Work Group. KDIGO clinical practice guideline for the diagnosis, evaluation, prevention, and treatment of chronic kidney disease-mineral and bone disorder (CKD-MBD). *Kidney Int.* 2009;76 (Suppl 113): S1-S130.
5. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

Immunodiagnostic Systems Limited,

10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK

Tel.: +44 191 519 0660

• Fax: +44 191 519 0760

e-mail: info.uk@idsplc.com

• www.idsplc.com

Immunodiagnostic Systems

USA Immunodiagnostic Systems (IDS) Inc,
948 Clopper Road, Gaithersburg, MD 20878, USA.

Tel.: +1 877 852 6210

• Fax: +1 301 990 4236

e-mail: info.us@idsplc.com

Franciaország Immunodiagnostic Systems SAS,
153 Avenue D'Italie, 75013 Paris, France

Tel.: +33 140 770 450

• Fax: +33 140 770 455

e-mail: info.fr@idsplc.com

Belgium Immunodiagnostic Systems S.A.,
Rue Ernest Solvay 101, 4000 Liège, Belgium

Tel.: +32 425 226 36

• Fax: +32 425 251 96

e-mail: info.be@idsplc.com

Németország Immunodiagnostic Systems GmbH (IDS GmbH)
Mainzer Landstrasse 49, 60329 Frankfurt am Main.

Tel.: +49 693 085 5025

• Fax: +49 693 085 5125

e-mail: info.de@idsplc.com

Skandinávia Immunodiagnostic Systems Nordic a/s (IDS Nordic a/s)
International House, Center Boulevard 5,

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EUA Immunodiagnostic Systems (IDS) Inc,
França Immunodiagnostic Systems SAS,
Bélgica Immunodiagnostic Systems S.A.,
Alemanha Immunodiagnostic Systems GmbH (IDS GmbH),