REF IS- 2030

IVD

1. Intended Use

For *In Vitro* Diagnostic Use

The IDS-iSYS 1,25 VitD^{xp} Control Set is for *in vitro* diagnostic use, for the quality control of the IDS-iSYS 1,25 VitD^{xp} on the IDS-iSYS Multi-Discipline Automated System.

Rx Only

2. Summary and Explanation

The regular use of control samples at several analyte levels is advised to ensure day-to-day validity of results. The IDS-iSYS 1,25 VitD^{xp} Control Set contains buffered human serum at two concentrations of 1,25(OH)₂D. The controls are used for monitoring the accuracy and precision of the IDS-iSYS 1,25 VitD^{xp} assay.

3. Warnings and Precautions

The IDS-iSYS 1,25 VitD^{xp} Control Set is for *in vitro* diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in the Instructions for Use. Immunodiagnostic Systems Limited (IDS) will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of non-compliance with the instructions provided.

CAUTION: This kit contains material of animal origin. Handle kit reagents as if capable of transmitting an infectious agent. Appropriate precautions and good laboratory practice must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV I and II), Hepatitis B surface antigen, antibody to Hepatitis C, and found negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled according to Biosafety Level 2.

Sodium Azide

Controls contain sodium azide (NaN₃) >0.1% (w/w) (<1%).

Classification under CLP:

Acute Tox. 4

Aquatic Chronic 3

Hazard statements:

EUH032 Contact with acids liberates very toxic gas
 H302 Harmful if swallowed
 H412 Harmful to aquatic life with long lasting effects

Precautionary statements:

P264 Wash hands thoroughly after handling.
 P270 Do not eat, drink or smoke when using this product.
 P273 Avoid release to the environment.
 P301/310 IF SWALLOWED: Immediately call a POISON CENTER or doctor.
 P330 Rinse mouth.
 P501 Dispose of contents/container to hazardous or special waste collection point.

4. Handling Precautions

The IDS-iSYS 1,25 VitD^{xp} controls are lyophilised. Ensure the controls equilibrate to room temperature (18 - 25°C). Reconstitute immediately before use. Add 1.2 mL of distilled or deionised water to each bottle. Replace the stopper. Leave for 10 minutes to reconstitute. Vortex for 2-3 seconds; invert gently to ensure complete reconstitution. Care should be taken to avoid the formation of foam.

If Controls are to be used more than once, they must be frozen (-20 °C) in the original vial within 15 minutes of reconstitution. When re-using frozen Controls; thaw at room temperature, vortex for 2-3 seconds and use within 15 minutes. Aliquots should not be re-frozen.

5. Shelf Life and Storage of Reagents

Store the controls in an upright position at 2 to 8°C prior to use. Self-defrosting freezers are not recommended for reconstituted controls storage.

Reagents shelf life	Controls
Before opening at 2 - 8 °C	To the expiry date
After reconstitution at 2-8°C	6 hours
After reconstitution at -20°C or lower in the original vial	28 days
Freeze/thaw cycle(s)	2
On board the IDS-iSYS *	4 hours

* Continuous on board stability

6. Materials

Materials Provided

CTL1	Human serum buffer matrix containing 1,25(OH) ₂ D and sodium azide as preservative (<1.0%, w/w). 6 each of 2 concentration level, 1.2 mL.
CTL2	

Materials Required But Not Provided

System

IS-310400	IDS-iSYS Multi-Discipline Automated System
IS-2000	IDS-iSYS 1,25 VitD ^{Xp}
IS-CC100	IDS-iSYS Cuvettes
IS-CS100	IDS-iSYS System Liquid (Syst I)
IS-CW100	IDS-iSYS Wash Solution (Wash S)
IS-CT100	IDS-iSYS Trigger Set
IS-6010	IDS-iSYS Cartridge Check System (CCS)
IS-CSC105	Sample Cups (500 µL)
IS-DW225	Disposable waste bags
IS-DS100	D-SORB solution

Optional equipment / materials

IS-2035	IDS-iSYS 1,25 VitD ^{Xp} Calibration Verifiers
IS-10DA	IDS-iSYS Diluent A

Disposable polypropylene 2 mL, 10.8 mm diameter conical skirted base, screw cap micro tubes and screw cap with O-ring (Sarstedt 72.609 and 65.716 or equivalent)

Precision Pipetting Devices

7. Assay

When using the Controls; pipette approximately 400 µL (218 µL per replicate) of each control into disposable 2 mL screw cap micro tubes and place on the System. Proceed according to the instructions of the IDS-iSYS Multi-Discipline Automated System User Manuals.

Discard the material in the micro tubes after use. DO NOT return material to the control vial.

All data required for the controls can be found on the mini CD. If the data for the lot of controls is not available on board the System, load the data using the mini CD provided with the control.

Use the IDS-iSYS 1,25 VitD^{Xp} Control Set for quality control. Controls should be tested minimally at least once every 24 hours when the test is in use and during every calibration in conformance with local, state and/or federal regulations or accreditation requirements and your laboratory's quality procedure. It is recommended that labs testing the specimens in multiple shifts in a day should measure the controls during each shift.

During a calibration, all control levels should be assayed in duplicate. The control values must be within the acceptable ranges specified. If a control is out of its specified range, the associated test results are invalid and samples must be retested. Recalibration may be required.

8. Symbols used

Catalogue Number

In Vitro Diagnostic Device

Manufacturer

Applied in accordance with directive 98/79/EC

Rx Only

Caution: US Federal Law restricts this device to sale by or on the order of a (licensed healthcare practitioner)

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660 • Fax: +44 191 519 0760
e-mail: info.uk@idsplc.com • www.idsplc.com

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Aquatic Chronic 3

Hazard statements:

EUH032 Contact with acids liberates very toxic gas

H302 Harmful if swallowed

H412 Harmful to aquatic life with long lasting effects

Precautionary statements:

P264 Wash hands thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P273 Avoid release to the environment.

P301/310 IF SWALLOWED: Immediately call a POISON CENTER or doctor.

P330 Rinse mouth.

P501 Dispose of contents/container to hazardous or special waste collection point.

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Aquatic Chronic 3

Hazard statements:

EUH032 Contact with acids liberates very toxic gas
H302 Harmful if swallowed
H412 Harmful to aquatic life with long lasting effects

Precautionary statements:

P264 Wash hands thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P273 Avoid release to the environment.
P301/310 IF SWALLOWED: Immediately call a POISON CENTER or doctor.
P330 Rinse mouth.
P501 Dispose of contents/container to hazardous or special waste collection point.

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EUA Immunodiagnostic Systems (IDS) Inc,
França Immunodiagnostic Systems SAS,
Bélgica Immunodiagnostic Systems S.A.,
Alemanha Immunodiagnostic Systems GmbH (IDS GmbH)