Seronorm™

Trace Elements Whole Blood L-1

REF 210105

LOT 1702821

2022-04

10 x 3 mL

IVD

Collaborating laboratories

The analytical data of Seronorm™ Trace Elements Whole Blood have been elaborated in collaboration with the following independent laboratories:

- ALS Scandinavia AB, Luleå, Sweden

EN

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C. No claims for MeHg.

After reconstitution, all elements are stable for 7 days when stored tightly capped at 2-8 °C.

Vials aliquoted and frozen within 30 minutes at ≤ -20 °C in screw cap tubes with O-ring are stable for 1 month.
Discard the remaining material after each use.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

Sample preparation before analysis

Samples were prepared by closed vessel MW-assisted digestion with concentrated nitric acid.

For MeHg measurement, please refer to Baxter, Rodushkin et al «Clinical Chemistry 53:1 000-000 (2007)»

Assignment of values

Table 1: The mean analytical values are derived from replicate analyses obtained through collaboration with independent laboratories and they are specific for this lot. The uncertainty of the analytical values is presented as a single number U (with a coverage factor $k = 2$), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given element. The analytical uncertainty (95 % confidence interval) is calculated from the analytical value \pm U. The analytical value for each element may vary due to differences in sample preparation, instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range.

The value assignment procedures have been established in accordance with the ISO 17511¹⁾ International standard using reference method procedures traceable to primary international standards.

Table 2: The approximate values should only be considered as an indication of the analyte concentration. No precise analytical values have been established for the elements reported in table 2 and lot to lot variation for these elements may occur.

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1) ISO 17511. In Vitro Diagnostic Medical Devices – Measurement of quantities in samples of biological origin – Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.

Explanation of abbreviations

ICP-SFMS: Inductively Coupled Plasma-Sector Field Mass Spectrometry

ICP-AES: Inductively Coupled Plasma-Atomic Emission Spectrometry

GC-ICP-SFMS: Gas Chromatography-Inductively Coupled Plasma-Sector Field Mass Spectrometry

Additional approximate values for Seronorm™ Trace Elements Whole Blood L-1 LOT 1702821

Element		Approximate value	Method
Barium	Ba	427 µg/L	ICP-SFMS
Bromine	Br	753 µg/L	ICP-SFMS
Cerium	Ce	75 ng/L	ICP-SFMS
Cesium	Cs	2,3 µg/L	ICP-SFMS
Dysprosium	Dy	10 ng/L	ICP-SFMS
Erbium	Er	7 ng/L	ICP-SFMS
Europium	Eu	<5 ng/L	ICP-SFMS
Gadolinium	Gd	13 ng/L	ICP-SFMS
Gallium	Ga	48 ng/L	ICP-SFMS
Gold	Au	<2 ng/L	ICP-SFMS
Hafnium	Hf	5 ng/L	ICP-SFMS
Holmium	Ho	<2 ng/L	ICP-SFMS
Iridium	Ir	<0,5 ng/L	ICP-SFMS
Iron	Fe	357 mg/L	ICP-SFMS
Iron	Fe	334 mg/L	ICP-AES
Lanthanum	La	73 ng/L	ICP-SFMS
Lutetium	Lu	<1 ng/L	ICP-SFMS
Neodymium	Nd	66 ng/L	ICP-SFMS
Niobium	Nb	30 ng/L	ICP-SFMS
Phosphorus	P	207 mg/L	ICP-SFMS
Phosphorus	P	203 mg/L	ICP-AES
Platinum	Pt	3,8 ng/L	ICP-SFMS
Potassium	K	1153 mg/L	ICP-SFMS
Potassium	K	1089 mg/L	ICP-AES
Praseodymium	Pr	17 ng/L	ICP-SFMS
Rhenium	Re	1,4 ng/L	ICP-SFMS
Rubidium	Rb	1,42 mg/L	ICP-SFMS
Samarium	Sm	12 ng/L	ICP-SFMS
Sodium	Na	1679 mg/L	ICP-SFMS
Sodium	Na	1598 mg/L	ICP-AES
Sulfur	S	1032 mg/L	ICP-SFMS
Sulfur	S	957 mg/L	ICP-AES
Tantalum	Ta	<1 ng/L	ICP-SFMS
Tellurium	Te	7 ng/L	ICP-SFMS
Terbium	Tb	1,3 ng/L	ICP-SFMS
Thorium	Th	9 ng/L	ICP-SFMS
Thulium	Tm	<1 ng/L	ICP-SFMS
Uranium	U	0,13 µg/L	ICP-SFMS
Ytterbium	Yb	8 ng/L	ICP-SFMS
Yttrium	Y	65 ng/L	ICP-SFMS
Zirconium	Zr	0,41 µg/L	ICP-SFMS