

# **VALIDATE®**

# Linearity and Calibration Verification Products Available for INSTRUMENTATION LABORATORY ACL TOP®

C€

# **HEMOSTASIS**

Five level, liquid-stable products specifically designed to challenge the IL ACL TOP system's full reportable range and eliminate manual dilutions.



#### **VALIDATE® D-DIMER**

ORDER NO. 902il

C€ ORD

**Typical Recovery:** 

200 - 3,250 ng/mL DDU

Matrix:

Human Plasma

Form:

3.0 mL, Liquid-Frozen

Storage:

-10 to -25° C

Open-Vial:

4 freeze-thaws through expiration

**VALIDATE® HEPARIN** 

ORDER NO. 903il

**Typical Recovery:** 

0.08 - 1.8 IU/mL

Matrix:

Human Plasma

Form

3.0 mL, Liquid-Frozen

Storage:

-10 to -25° C

Open-Vial:

2 freeze-thaws through expiration

VALIDATE® FIBRINOGEN

ORDER NO. 904il

C€

**Typical Recovery:** 

50 - 850 mg/dL

Matrix:

Human Plasma

Form:

2.0 mL, Liquid-Frozen

Storage:

-10 to -25° C

Open-Vial:

Single-use, 6 hours

# EASY • FAST • EFFICIENT



## Liquid, Ready-to-Use

Increase productivity and avoid the time, error and inconvenience of reconstitution. Supplied in dropper bottles, simply bring to room temperature and dispense into the sample cup, reducing the risk of contamination.



# **Open-Vial Stability**

Freeze-thaw cycles allow for multiple uses cap the bottle, return to storage and use again. For Heparin and Fibrinogen, order once and receive a second, no charge shipment, in 6 months.



## **Linearity & Peer Groups**

5 Level products with free data reduction services offers instant linearity results via MSDRx with ready access to Peer Statistics.



# Accreditation

For use at installation, or during the laboratory's routine 6-month testing interval, rely on VALIDATE to satisfy ISO 15189, CAP, JCI and local regulatory requirements, including AMR validation.



#### Quality

Products are U.S. FDA cleared or listed and are CE marked. Our facility is registered with the U.S. FDA and our QMS is certified to ISO 13485.



## **Troubleshooting**

VALIDATE is ideal for use after preventative maintenance or to troubleshoot across the full reportable range.



mainestandards.com / 800-377-9684 or 207-892-1300 / CDx-CustomerService@lgcgroup.com

