

TECHNICAL SHEET IDYLLA™ IDH1-2 MUTATION ASSAY KIT



The **Idylla™ IDH1-2 Mutation Assay Kit** performed on the Biocartis Idylla™ System, is a fully automated in vitro diagnostic assay for the qualitative detection of five IDH1 mutations in **codon R132**, four IDH2 mutations in **codon R140** and six IDH2 mutations in **codon R172**.

The Idylla™ IDH1-2 Mutation Assay Kit uses extracted DNA from blood (EDTA) or human bone marrow (EDTA) from patients with Acute Myeloid Leukemia (AML) or DNA extracted from FFPE tissue from patients with glioma. The Kit is also compatible with FFPE samples from patients with glioma. The Assay is a **sample-to-result** real-time PCR, intended for research use only.

FEATURES

IDH1 Mutation Detection

Codon 132	R132C	c.394C>T
	R132H	c.395G>A
	R132G	c.394C>G
	R132S	c.394C>A
	R132L	c.395G>T

IDH2 Mutation Detection

Codon 140	R140Q	c.419G>A
	R140L	c.419G>T
	R140G	c.418C>G
	R140W	c.418C>T
Codon 172	R172K	c.515G>A
	R172M	c.515G>T
	R172G	c.514A>G
	R172S	c.516G>T
	R172S	c.516G>C
	R172W	c.514A>T

Specimen requirements

Sample Type	Extracted DNA: 50 µl (concentration 10 µl/g) OR FFPE: It is recommended to use maximum 3 FFPE tissue sections.
Neoplastic cells	Extracted DNA: DNA should be extracted from samples with ≥ 10% neoplastic cells. If less, an enrichment procedure is required. OR FFPE: ≥ 10% neoplastic cells. If less, macro-dissection is required.

Supported blood and bone marrow collection tubes

EDTA tubes

Total turnaround time

Time 95 minutes

Result reporting

Report Qualitative genotype call + Cq Value

Catalog numbers

Idylla™ IDH1-2 Mutation Assay (Vial) A0181/6

Idylla™ DNA Cartridge A0191/6

The Idylla™ IDH1-2 Mutation Assay Kit consists of a Cartridge and a Vial.



Idylla™ IDH1-2 Mutation Assay Kit is for Research Use Only (RUO), not for use in diagnostic procedures. Idylla™ Platform is available as an IVD in Europe, the US and many other countries. Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis.

Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with the local Biocartis representative.
© July 2023, Biocartis NV. All rights reserved.

WWW.FLO56.EN.N.07/2023

Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx