

Epredia Cytospin4™

Adopt **Efficiency**. Maintain **Integrity**.

Trusted tools that **enhance precision cancer diagnostics**.

Cytology solutions from Epredia have been the standard in liquid based cytology for over 50 years.

Maintain cellular integrity and carefully separate and deposit cells in a monolayer. The Cytospin system is adaptive and convenient for use in cytology, virology, microbiology, haematology, oncology and research laboratories.



epredia
Enhancing precision
cancer diagnostics

The ultimate thin-layer cell preparation system

Cyto centrifugation provides an efficient way to present cells from a liquid matrix for microscopic examination - maximising cell yield while maintaining morphological integrity.

Optimized Cell Yield

The Cytospin offers a quick, efficient and cost-effective technique to provide superb cell recovery from hypocellular samples with advanced preservation of cell morphology, while preserving precious sample for further testing or preparation of cell blocks.

Excellent Nuclear Presentation

Using centrifugal principles, the Cytospin deposits a monolayer of cells onto a clearly defined area of a glass slide. Constructively flattening cell allows for excellent nuclear presentation during subsequent microscopic analysis or digital scanning, with a good representation of all cell types present within the homogenous liquid sample in a uniform layer.

Safety and Simplicity

Designed to protect the operator before, during and after sample preparation, the Cytospin features an autoclavable sealed head which can be loaded and unloaded in a microbiological safety cabinet to reduce the chance of exposure to any potential pathogens. The Cytospin offers an extremely simple technique to concentrate up to 12 samples in one run. The intuitive control panel provides a simple interface for inputting time, speed and acceleration settings, with a protected program memory for up to 23 commonly used protocols.

Cytospin4 and Accessories

Cytospin4 Cyto centrifuge Unit Only	Catalogue No.
-------------------------------------	---------------

100-120V, US Power Cord	A78300003
-------------------------	-----------

220-240V, UK Mainslead Set	A78300002
----------------------------	-----------

220-240V, Schuko Mainslead Set	A78300004
--------------------------------	-----------

220-240V, Chinese Mainslead Set	A78300005
---------------------------------	-----------



Enhance the use of your Cytospin centrifuge with a wide variety of high quality EpreDia consumables, including funnels, slides and fixation solutions.

Contact your EpreDia Representative today for more information, or visit our website at www.epredia.com/cytology-solutions


LiveLab

Experience the Cytospin for yourself, safely and conveniently with the **EpreDia LiveLab**. If you are interested in a LiveLab session, please reach out to your EpreDia Representative.


ePredia

Enhancing precision
cancer diagnostics

© 2021 EpreDia Group. All rights reserved. All trademarks are the property of EpreDia and its subsidiaries unless otherwise specified. M51016 R092021

Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

