

Helicobacter pylori Detection



- ***CoproELISA™ H. pylori***
- ***CoproStrip™ H. pylori***



***H. pylori* Stool-Antigen Immunoassays to Fit Your Testing Requirements**

Savyon Diagnostics is now offering two assays for detection of *H. pylori* active infection. The CoproELISA™ *H. pylori* test which allows manual or automatic high-throughput testing within 75 min and the CoproStrip™ *H. pylori* rapid test which offers simple, one-step testing within 10 min.

Both tests are based on the detection of a unique *H. pylori* antigen from feces of symptomatic patients.

Current consensus guidelines* recommend testing *H. pylori* active infection by stool antigen immunoassays such as the CoproStrip™ and CoproELISA™ *H. pylori*

Helicobacter pylori is now recognized as one of the most common and medically important pathogens worldwide. It has been firmly established as an etiologic agent in chronic gastritis and peptic ulcer disease, and has been associated with Mucosa-Associated Lymphoid Tissue Lymphoma (MALT) and Gastric Adenocarcinoma. The stool antigen test is now accepted as an accurate, non-invasive test to be used both pre- and post antimicrobial treatment.

* The American College of Gastroenterology (2005) and the American Gastroenterological Association (2007)

CoproELISA™ *H. pylori*

- ❖ Accurate (>95% sensitivity & specificity)
- ❖ Flexible (1-96 samples /run)
- ❖ Short & Simple Protocol
- ❖ Manual or Automation Procedures
- ❖ Affordable



**Easy to use,
Easy to read**

CoproStrip™ *H. pylori*

- ❖ Rapid (10 min)
- ❖ Accurate (>95% sensitivity & specificity)
- ❖ User friendly
- ❖ Only 1 min of set-up
- ❖ Affordable



Test Name	Tests/kit	Catalog No.
CoproELISA™ <i>H. pylori</i>	96	774-01
CoproStrip™ <i>H. pylori</i>	20	41221

3542 01-05/14