



AccuSet™ SARS-CoV-2 Performance Panel

ACCELERATE SEROLOGY ASSAY DEVELOPMENT, VALIDATION AND IMPLEMENTATION

ABOUT SERACARE

TRUSTED SUPPLIER TO DIAGNOSTIC MANUFACTURERS AND CLINICAL LABS FOR OVER 30 YEARS

HIGH-QUALITY CONTROL PRODUCTS, RAW BIOLOGICAL MATERIALS, AND IMMUNOASSAY REAGENTS

INNOVATIVE TOOLS AND TECHNOLOGIES TO PROVIDE ASSURANCE IN DIAGNOSTIC ASSAY PERFORMANCE AND TEST RESULTS

FOR MORE INFORMATION, PLEASE VISIT OUR WEBSITE: WWW.SERACARE.COM

AccuSet™ SARS-CoV-2 Performance Panel is intended to provide an out-of-the-box solution to evaluate SARS-CoV-2 antibody detection assays with highly characterized human specimens. The eleven-member validation panel consists of undiluted, naturally occurring human-plasma samples. Ten panel members represent bleeds from multiple individuals positive for antibodies to SARS-CoV-2; a single negative member is also included. Test results from commercially available SARS-CoV-2 antibody assays are included for comparative analysis.

Product Description	Material Number	Configuration	Fill Size
AccuSet SARS-CoV-2 Performance Panel	0820-0410	10 positive members 1 negative member	1.0 mL

Not for In Vitro Diagnostic Use. Research Use Only.

AccuSet SARS-CoV-2 Serology Solution

- **True Patient Specimens:** Unique, individual plasma samples confirmed positive for SARS-CoV-2 antibodies
- **Highly Characterized:** Comprehensive data sheet with results from leading commercially available assays for SARS-CoV-2 antibody detection
- **Out-of-the-box Solution:** Ideal for quickly and consistently performing assay validation, verification, and installation

Whether generating validation data for a regulatory submission or performing assay verification in a clinical laboratory setting, LGC SeraCare's panel products offer a comprehensive solution to rigorously evaluate assay performance. Our suite of products will allow you to confidently develop, challenge, validate, verify and troubleshoot your assays from initial proof of concept, to production, and into the clinical laboratory.

- **Save time and money** spent sourcing hard-to-find, uncharacterized specimens
- **Accelerate assay development, installation and verification** with pre-characterized, unique donor samples
- **Ensure test performance** with human disease state material

To learn more about the AccuSet SARS-CoV-2 Performance Panel and our complete line of SARS-CoV-2 Quality Solutions, visit www.seracare.com/SARSCoV2 or contact customer service at 800.676.1881

ORDERING INFORMATION

To place an order, please contact us at +1.508.244.6400 and 800.676.1881 or email customerservice@seracare.com.