



Chromascopics®

Urinalysis Control with Microscopics

- Specifically designed for Siemens Clinitek Novus, Atellica UAS 800 & Atellica 1500 AUS
- Microscopics sediment elements include: RBC, WBC, calcium oxalate dihydrate crystals and now introducing E. coli bacteria
- 24-month refrigerated shelf life from date of manufacture*



NEW CHROMASCOPICS Urinalysis Control with Microscopics may be used for validating the processing of urine samples prior to the microscopic evaluation of urine sediment. The control may also be used for confirmatory tests such as: K-Check, Ictotest and for β hCG screening methods. CHROMASCOPICS Urinalysis Control is designed for use with major urinalysis analyzers including: Clinitek Novus, Clinitek Advantus, Clinitek Atlas, Clinitek Status, Clinitek Status Plus, Clinitek 50, Clinitek 500, Atellica UAS 800 & Atellica 1500 AUS. CHROMASCOPICS is intended for use in traditional testing environments including reference labs, hospitals, doctors offices and clinical lab settings.

Our thoughtfully designed packaging includes barcodes on the bottles for lot numbers and expiration dates and we've made it easy to identify each level.

Dropper Analytes*

Bacteria, Bilirubin, Blood, Clarity, Color, Creatinine, Calcium Oxalate Crystals, Glucose, hCG, Ketones, Leukocytes, Microalbumin, Nitrite, pH, Protein, Red Blood Cells, Specific Gravity, Urobilinogen and White Blood Cells

Catalog #	Contents	Size
1540-01	Chromascopics (INT) – 2 Level Set	4 x 120 mL
1540-02	Chromascopics (US) – 2 Level Set	4 x 120 mL

Coming Soon: (310) 536-0006

Visit us at: quantimetrix.com

*Refer to the package insert of currently available lots for specific analyte, stability claims and instrument compatibility.

10-0169_0518 © 2018 Quantimetrix. All Rights Reserved.



Distributed by Abacus dx

1800 ABACUS (AUS) 0800 222 170 (NZ) | info@abacusdx.com | www.abacusdx.com

abacus dx