

Blood Collection Tube

Cyto-Chex® BCT is a direct-draw blood collection tube used for immunophenotyping of white blood cells by flow cytometry.



Features and Benefits

- Direct-draw blood collection tube for the preservation of whole blood specimens
- Preservative minimizes the adverse effects of time, storage and transport conditions on sample integrity
- Samples can be held for up to 14 days prior to analysis, allowing for batch processing

Cell Stabilization

- Maintains cellular morphology and surface antigen expression
- Samples are stable at room temperature for convenient transport and storage
- FDA 510(k) cleared for consistent recovery of HIV-associated lymphocyte subsets for 14 days

Batching

Samples collected in Cyto-Chex BCT can be stored at room temperature for 14 days prior to analysis, improving laboratory efficiency.

- Reduce valuable labor time
- Reduce daily flow cytometry start-up costs
- Reduce the need for patient redraws

Associated Product

Streck Cell Preservative™ is a liquid preservative that maintains cellular antigen expression, including cluster of differentiation (CD) markers and cell morphology of biological samples, for analysis by flow cytometry.

Ordering Information

Product Description	Catalog #
Cyto-Chex BCT Glass 5.0 mL	
6-tube pack	213361
25-tube box	213362
100-tube box	213386
Cyto-Chex BCT Glass 2.0 mL	
6-tube pack	213559
100-tube box	218980